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GFP originally comes from the jellyfish *Aequorea victoria*.

238 amino acids. Molecular weight 27 kDa. GFP fluoresces green when exposed to U.V. light (excitation 395 nm, emission 475 nm). The fluorophore is internal and consists of cyclisation of the sequence Phe₆₄-Ser-Tyr-Gly-Val-Gn₆₉. GFP exposes many hydrophobic amino acids residues.

The fluorescence properties of GFP have been changed by genetic engineering leading to several mutants, especially EGFP (Enhanced GFP, fluorophore = Leu₆₄-Thr-Tyr-Gly-Val-Gn₆₉) which have red-shifted excitation spectra (maximal excitation peak at 488 nm) and fluoresces (at 509 nm) 35 more brightly than wild-type GFP. EGFP extinction coefficient at 488 nm = 55 000 L.mol⁻¹.cm⁻¹

pET15-EGFP = MGSSHHHHHHSSGLVPRGSH MVSKGEELFTGVVPIVELDGDVNGHKFSVSGEGEGDATY GKLTLKFICT TGKLPVPWPTLVTTLTLYGVQCFSTRYPDHMKQHDFFKSAMP EGYVQERTIFKDDGNYKTRAEVKFEGDTLVNRIELKGIDFKEDGNILGH KLEYNYNSHNVYIMADKQKNGIKVNFKIRHNIEDGSVQLADHYQQNTPIGDGPVLLPDNH YLSTQSALSKDPNEKRDMV LLEFVTAAGI TLGMDELYK ; 20+239 = 259 aa : **29,11 kda** ; **pi calculé = 5,83**.

LVPRIGS = THROMBIN CLEAVAGE SITE SEQUENCE

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Annexe 7. Donn es techniques concernant le support His Trap FF

His Trap FF, 1 mL column , TECHNICAL SPECIFICATIONS

Medium : Ni Sepharose 6 Fast Flow

Column volume : 1 ml

Column dimensions : 0.7 x 2.5 cm

Dynamic binding capacity: Approx. 40 mg histidine-tagged protein/ml medium (*Protein binding capacity is protein-to-protein dependent*)

Recommended flow rate : 1 ml/min

Max. pressure 3 bar (0.3 MPa, 42 psi)

pH stability and Compatibility : Stable in all commonly used buffers, reducing agents, denaturants, and detergents

Storage 20% ethanol Storage temperature 4 C to 30 C

Notice : 11-0008-88AE

[http://www5.gelifesciences.com/aptrix/upp00919.nsf/Content/62743F04C34D5EBBC1256F24000846BD/\\$file/11000888AE.pdf](http://www5.gelifesciences.com/aptrix/upp00919.nsf/Content/62743F04C34D5EBBC1256F24000846BD/$file/11000888AE.pdf)

<p>His trap FF Stripping and recharging</p> <p>Note: The column does not have to be stripped and recharged between each purification if the same protein is going to be purified; it is sufficient to strip and recharge it after approximately 5 purifications, depending on the cell extract, extract volume, target protein, etc.</p> <p>Recommended stripping buffer: 20 mM sodium phosphate, 0.5 M NaCl, 50 mM EDTA, pH 7.4</p> <p>Strip the column by washing with at least 5–10 column volumes of stripping buffer. Wash with at least 5–10 column volumes of binding buffer and 5–10 column volumes of distilled water before recharging the column.</p> <p>Recharge the water-washed column by loading 0.5 ml or 2.5 ml of 0.1 M NiSO₄ in distilled water on HisTrap FF 1 ml and 5 ml column, respectively. Salts of other metals, chlorides, or sulfates, may also be used (see Optimization). Wash with 5 column volumes distilled water, and 5 column volumes binding buffer (to adjust pH) before storage in 20 % ethanol.</p>	<p>His trap FF cleaning</p> <p>When an increase in backpressure is seen, the column should be cleaned. Before cleaning, strip off Ni²⁺ ions using the recommended procedure described above.</p> <p>After cleaning, store in 20 % ethanol (wash with 5 column volumes) or recharge with Ni²⁺ prior to storage in ethanol.</p> <p>The Ni²⁺-stripped column can be cleaned by the following methods;</p> <ul style="list-style-type: none"> Remove ionically bound proteins by washing with several column volumes of 1.5 M NaCl; then wash with approx. 10 column volumes of distilled water. Remove precipitated proteins, hydrophobically bound proteins, and lipoproteins by washing the column with 1 M NaOH, contact time usually 1–2 hours (12 hours or more for endotoxin removal). Then wash with approx. 10 column volumes of binding buffer, followed by 5–10 column volumes of distilled water. Remove hydrophobically bound proteins, lipoproteins, and lipids by washing with 5–10 column volumes 30 % isopropanol for about 15–20 minutes. Then wash with approx. 10 column volumes of distilled water. <p>Alternatively, wash with 2 column volumes of detergent in a basic or acidic solution. Use, for example, 0.1–0.5% nonionic detergent in 0.1 M acetic acid, contact time 1–2 hours. After treatment, always remove residual detergent by washing with at least 5 column volumes of 70 % ethanol. Then wash with approx. 10 column volumes of distilled water.</p>
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------